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The Buccal Activator is characterised by a double channel that cover the full dental arches, a vestibular band, lingual and an occlusal base that create a correct dental arch (angle class I).

Principle

The active mouth orthosis acts as a dental guide. This dynamic, so-called functional device is a perfection proprioceptive simulation (teeth, ligaments, muscles,...), and supports balanced muscular work and, thus, huge relief for the temporomandibular joints (jaw).

Goals

Through natural muscular power and perfect proprioceptive simulation, this functional guide enables:
- an improvement to the chewing movement,
- mandibular repositioning,
- balanced occlusion,
- real effectiveness against bruxism and temporomandibular problems,
- re-education of chewing, breathing through the nose, swallowing and phonation.
- reduced snoring and better quality sleep.

Remark: the Buccal Activator can also be used in traditional orthodontics before a permanent device is fitted, or as a retainer after treatment.

Functions

The Buccal Activator has a double active function:
- Because of its elastic material, with shape-retention memory, it works like an active spring, forcing the muscles to return to their primary functions.
- Because of conscious chewing, patients must work actively with the Buccal Activator (self-stimulation).

Postural action

The Buccal Activator improves the body's posture due to its descending chain action. First, this dynamic orthosis is designed to rebalance the temporomandibular joints (jaw). Just like the feet and the eyes, the jaw is a strategic part of the human body. If it works poorly, over time, it can cause descending chain postural problems. It is common to find migraines, cervical pain, back pain and even lumbar pain associated with jaw problems.

Use

The Buccal Activator can be used during the day to do maxillofacial strengthening exercises that guide the jaw into the right position, and at night to ensure a perfect loosening of the temporomandibular joints (jaw). In particular, the orthosis is used in orthodontics. It is surprising to note, that by actively using this orthosis, the teeth find a balance in the mouth by placing themselves correctly vis-à-vis each other.

How to use your Buccal Activator (for children and adults) ?

The Buccal Activator is a functional, self-stimulation device. First, it must be worn repetitively during the day for seven to ten days. In the case of nocturnal bruxism (clamping the teeth at night), after this period of daily wear it is recommended to wear them during the day (intermittent active use) and at night (continuously, resting). As regards the daily exercises, you should try to keep your lips together and not bite down on the device but chew it (biting hard will not improve results, quite the contrary!). Furthermore, you should not speak with the Buccal Activator in your mouth.

Keep to the adaptation time for wearing the Buccal Activator. It may be uncomfortable for the first few days. That is why it is recommended you go slow, starting with a few minutes per day, building up to between twenty and thirty minutes (it varies from case to case and only a professional can advise you different times).

The Buccal Activator should not cause pain or injuries. If it does, ask your dentist to adjust the device to reduce rubbing against the mucous membrane. Optimum oral hygiene is required when wearing the Buccal Activator. If you feel pain (either in the joints or muscles) at the beginning of the treatment, it is best to stop as soon as the pain arises and to start the exercises again a little later (you must listen to your body and allow it to rest if there is a problem).

For those who need to wear the device at night: do not worry if you do not manage to wear it all night at first, especially if you breathe through your mouth. Repeating the exercises regularly is key. The incisors may be sensitive at the start of treatment, but this will quickly disappear. Hypersalivation is entirely normal when you start wearing the Buccal Activator. This will reduce gradually and you will be able to wear it continuously without problems.

If the top of the Buccal Activator is too big and irritates your gums, you can trim a couple of millimetres. The lateral branches can also be cut off in case of discomfort (or in the case they cause nausea).

How to achieve a Stimulating Postural Action (SPA) :

During the day, facilitate the work of the chewing muscles and decrease the constraints on the temporomandibular joints (jaw). Relaxing the jaw joint instantly relieves tension in the muscles of the neck and nape of the neck. The head will become more balanced and the constraints caused by an unbalanced jaw will become less marked in descending chain. Wear your Buccal Activator when walking or playing sport. It is recommended you chew the orthosis to ensure the jaw relaxes more. For intense physical efforts where you must breathe through the mouth, it is recommended you only wear the orthosis during warm-up. To achieve a proprioceptive balance effect, it is recommended you wear the Buccal Activator for 5 to 10 minutes per day. For greater muscle relaxation, it is recommended you perform a couple of stretching exercises after using the Buccal Activator for a few minutes.

1. Turn the head left then right and bow it on each side.
 2. Turn your torso (including the shoulders) to each side and bow it on each side.
 3. Keeping your legs straight, lean forward as if you wanted to touch the floor with your hands.
 4. Crouch by bending at the knees then straighten back up.
- Repeat the exercise two or three times as you wish. Your tensions should be relieved immediately.

The Buccal Activator should be :

- Rinsed with lukewarm water after each use
- Brushed with toothpaste once per day
- When not in use, it should be stored in its box with the lid open.

What is the Buccal Activator's lifespan ?

Because this therapeutic tool is made of silicone, it will lose its elasticity and performance over time. In general, the Buccal Activator's lifespan is between 3 and 6 months. It depends on a number of factors such as level of activity, chewing force and care of the orthosis.

Reasons not to wear the Buccal Activator ? Significant tooth movement is a reason not to wear the Buccal Activator. If you can feel your teeth moving by pressing them with your fingers, do not chew the Buccal Activator. If in doubt, ask your dentist for advice. If you have just completed orthodontic treatment it is advised to wait 3 months before using the Buccal Activator and to ask your dentist for advice. If you have a bridge, implant or crowns, it is recommended you chew the Buccal Activator no more than 5 to 10 minutes per day. If your goal is to improve your teeth's position, it is best to be monitored by a dentist specialised in functional orthodontics.

Wearing the orthosis is contraindicated :

- if one of your incisor is clearly longer than others (< 3mm)
- in the case of pain in the jaws, it is advised to wear the orthosis passively only (without any exercises, resting) until the pain disappear.

Material : silicon without bisphenol nor phthalates

All of the information regarding the Buccal Activator, whether on our website or in the catalogue, is for information only. Neither it nor the Tedop products are sold or designed to diagnose, treat or heal any problem. Contact your healthcare professional for more detailed information, diagnosis, screening, treatment and advice on how to use this information or Tedop products.

Precaution for use

It is recommended to clean the orthosis before use. For people with import mouth misformations (retrognathia our advanced prognathism) or oral pathologies, please consult your doctor to eventually get the dental range of products.

If you are currently under treatment for implants, prosthetic or orthodontic, please ask your dentist advice before using the orthosis.

If you are affected by periodontitis (teeth loosening), ask advice to your dentist. This orthosis should not be used by people suffering from epilepsy or convulsions.

Side effects

In the case of irritation of the palate and hypersensitivity of the teeth, it is recommended to reduce the duration of use. In the case of a lesion of the palate, it is preferable to stop wearing the orthosis and ask advice to your doctor.

In the case of nausea and vomiting, it is possible that the posterior branches may be too long. You can make them shorter with scissors. If nausea and vomiting still occur, it is possible that your palate is hypersensitive. Please inform your dentist for a desensitization treatment.

Slight gum bleeding can occur in the beginning when you wear the orthosis. Progressive wearing is advised to ensure an reinforcement of the tissues and a stop to the bleedings.

The orthosis should in no case cause wounds or pains. If this occur, ask your practitioner to modify the orthosis to reduce the friction with the mucous membrane. Wearing the orthosis require an optimal oral hygiene. In case of pain (muscular or articular) in the beginning of your treatment, it is preferable to stop using the orthosis as soon as pain appear et start exercises again later (It is essential to listen to your body and let it rest should pain occur)

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